

# EasyTest Patch Kit

**EPK**

CS 1000  
CS 1939  
CSI-C-11  
HY-TRAX™  
RBSA  
CSM  
TFL  
TFH  
FCU  
MCS  
AS  
SMU  
CTU  
**EPK**

## Features and Benefits

- User-friendly visual analysis of solid contamination
- Compatible with mineral-based hydraulic fluids and lubricants, and petroleum distillates
- Enables solid contaminant quantification and identification
- Provides on-site results in a matter of minutes

## Applications

- Perform quick on-site determination of contamination levels of solid particles
- Supplement on-site laboratories
- Use as a tool to demonstrate need for solid contamination mitigation



The Schroeder EasyTest Patch Kit (EPK) provides the necessary tools to determine the level of solid particulate contamination present in a fluid sample. Using the vacuum pump contained in the kit, the fluid sample is drawn through a membrane patch. The residual particulate left on the patch is viewed under a microscope and compared to photos of known contamination levels in the L-2711 Schroeder Contamination Handbook (included).



## Description

Trouble  
Check Plus  
HMG2500  
HMG4000



Schroeder:  
EasyTest Patch Kit (EPK)

7630322



## Model Selection + Items Supplied P/N 7640674

HTB  
RFSA  
HFS-BC  
HFS-15  
MFD-BC  
MFS, MFD  
HY-TRAX™  
Retrofit System  
MFD-MV  
MFS-HV  
AMS, AMD  
FS  
AMFS

Kit as supplied includes:

Quantity	Description	P/N
1	Hand-held vacuum pump	7619502
3	Syringe, 30 mL	7626475
50	Disposable Petri Dishes	7630320
1	Forceps	7626481
1	Membrane patches, 0.45 µm, 25 mm, (100 pack)	2701997
1	Membrane patches, 0.8 µm, 25 mm, (100 pack)	2701952
1	Carrying Case	7640195
1	Microscope, 10x - 200x	7635242
1	Plastic funnel	7626479
1	Solvent dispenser bottle	7626473
1	Solvent dispenser bottle cap	7640496
3	Plastic sample bottle, 4 oz.	7626480
1	Solvent patch holder	7632471
1	Tubing, Tygon 3"	7624738
1	10' section of ¼" LDPE tubing	2701999
1	L-2711 Contamination Handbook & Instructions	7627179

### NOTES:

Solvent is not supplied w/ the EPK. Recommended solvents include Heptane (99% by GLC), or Isopropyl Alcohol.

Kit contents are subject to change at the discretion of the manufacturer.

LSN, LSA, LSW  
X Series  
OLF Compact  
OLF  
OLF-P  
NxTM  
VEU  
IXU  
Triton-A  
Triton-E  
NAV  
SVD01  
SVD